

|               |   |
|---------------|---|
| Product Name: | LRP-HC-72-TS  |
| Description:  | <p>TempLog Plus Series Glass Door Laboratory Refrigerator<br/> 72 Cu. Ft. capacity<br/> Three swing glass doors, self-closing (energy efficient, double pane vacuum insulated)<br/> Microprocessor temperature controller<br/> Adjustable operating temperature range: 1°C - 10°C<br/> Touch Screen Digital Temperature Display and Alarm Module with battery back-up to support temperature monitoring and alarms<br/> Password-protected control panel<br/> Sample simulated digital temperature display; .1° Resolution, adjustable temperature display in Celsius or Fahrenheit<br/> Active temperature scale bar display<br/> Resettable 1-point min/max history<br/> Digital 24-hour temperature chart<br/> Temperature data logging<br/> USB port for data transfer and software updates<br/> Real-time operational display clock<br/> 4-20mA output<br/> Electronic troubleshooting guide<br/> Electronic QR Code to access Owner's Manual &amp; technical support contact information<br/> Adjustable alarm volume<br/> Audible and visual high/low temperature alarms<br/> Sensor error alarm<br/> Power failure alarm<br/> Door ajar alarm<br/> Low battery alarm<br/> Operation indicator<br/> Alarm validation (audible, visual, and remote alarm)<br/> Alarm mute with ring back<br/> Remote alarm contacts<br/> Twelve adjustable shelves (adjustable in ½" increments)<br/> Forced draft air circulation<br/> Keyed door lock<br/> Magnetic door gasket for positive seal<br/> Heavy duty door handles<br/> LED interior lights are safety shielded and switch controlled<br/> Swivel casters (adds 4" to height)<br/> Probe access port (3/8")<br/> White laminated vinyl exterior<br/> White aluminum liner and stainless steel floor interior<br/> HFC-free refrigerant (R290), foam insulation and packaging<br/> Three year parts and labor warranty, plus an additional two year compressor parts warranty<br/> Exterior dimensions: 78 1/4"W x 31"D x 82 5/8" H with casters; 78 5/8" H without casters<br/> Cycle defrost<br/> 115V, 60 Hz, 8.5 Amps, 1/2 HP, BTU/HR 7315<br/> UL/C-UL listed<br/> ENERGY STAR CERTIFIED<br/> Shipping weight: 690 lbs.</p> |