

Product Name:	LRP-HC-RFC1040
Description:	<p>Plus Series Laboratory Refrigerator &amp; Freezer Combo (-40°C)  Freestanding  5.2 cu. ft. capacity (Refrigerator) / 4 cu. ft. capacity (Freezer)  One solid swing door (fridge), one solid swing door (freezer)  Microprocessor temperature controller  Digital temperature display  Audible and visual high and low temperature alarms  2 probes; 1 in air and 1 in sample bottle  Remote alarm contacts  Adjustable operating temperature range- Freezer: -37°C – -43°C; Refrigerator: 2°C to 8°C  3 adjustable shelves (freezer); 3 adjustable shelves (fridge)  Bottom footed shelf to allow air circulation  Back wall evaporator cover  Keyed door locks  Magnetic door gasket for positive seal  Recessed door handles  Stacking kit  Leveling legs  LED interior light  Probe access port (3/8")  White cabinet and door  High density urethane foam cabinet and door insulation  HFC-free refrigerant (R600a) Refrigerator; Natural, Hydrocarbon refrigerant (R290) Freezer  Two year parts and labor warranty, plus an additional three year compressor parts only warranty  Exterior dimensions: 23 3/4" W x 24"D x 65"H  <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i>  Manual defrost (freezer); Cycle defrost (refrigerator)  115V, 60 Hz  2.2 Amps, 1/4 HP (freezer); 1.3 Amps, 1/6 HP (refrigerator)  UL/C-UL listed  Shipping weight: 268 lbs.</p>