

Product Name:	LRP-HC-RFC9G
Description:	<p>Plus Series Refrigerator & Freezer Combination 9 Cu. Ft. Capacity (4.6 Ref./4.2 Freezer) Two exterior swing doors (glass door refrigerator/solid door freezer), right hinged* Microprocessor temperature controller, refrigerator and freezer Digital temperature display, refrigerator and freezer Adjustable operating temperature range: Refrigerator: 1°C – 10°C; Freezer: -15°C – -25°C Audible and visual high and low temperature alarms, refrigerator and freezer 2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer Remote alarm contacts, refrigerator and freezer 2 adjustable/1 fixed shelves refrigerator with guard rail on back to prevent product from being pushed against back wall/2 fixed shelves freezer Back wall evaporator cover, Refrigerator Forced draft circulation, Refrigerator Keyed door locks Magnetic door gasket for positive seal Recessed door handle, freezer; heavy duty door handle, refrigerator Leveling legs Stacking kit LED interior light, Refrigerator Two probe access ports (3/8"); one, refrigerator/one, freezer White powder coated exterior High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Exterior dimensions: 23 3/4" W x 25 1/2"D x 66 3/4"H Two year parts and labor warranty, plus an additional three year compressor parts warranty Manual defrost freezer / cycle defrost refrigerator 115V, 60 Hz, 1.74 Amps Refrigerator and 1.5 Amps Freezer, 1/8 HP Refrigerator and Freezer UL/C-UL listed Shipping weight: 282 lbs. *Please note this unit requires 2 outlets. *This product is available left hinged, product number LRP-HC-RFC9G-LH</p>