

Product Name:	LRP-HC-SLB-49
Description:	<p>Basic Series Solid Door Laboratory Refrigerator 49 Cu. Ft. capacity Two swing solid doors, self-closing, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in sample bottle Remote alarm contacts Adjustable operating temperature range: 1°C – 10°C Keyed door locks Magnetic door gasket for positive seal Eight adjustable shelves (adjustable in ½” increments) Forced draft air circulation Heavy duty door handles LED interior lights are safety shielded and switch controlled Swivel Casters (3 ½”) probe access port (⅜”) White powder-coated steel exterior and interior Hydrocarbon, natural refrigerant (R290) One year parts and labor warranty, plus an additional four year compressor parts warranty Exterior dimensions: 54” W x 34 ¾” D x 81 ¾” H with casters; 78 ¼” H without casters Cycle defrost 115V, 60 Hz, 4.5 Amps, 1/4 HP ENERGY STAR CERTIFIED ETL listed Shipping weight: 446 lbs.</p>