

Product Name:	LRP-HC-UCBI-0420SS-LH
Description:	<p>Plus Series Undercounter Freezer Built-In-LEFT HINGED 4.2 Cu. Ft. Capacity One stainless steel swing door, left hinged Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in sample bottle Remote alarm contacts Adjustable operating temperature range: -15°C – -25°C Two fixed shelves Keyed door lock Recessed door handle Magnetic door gasket for positive seal Leveling legs Probe access port (3/8") Stainless steel exterior High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Two year parts and labor warranty, plus an additional three year compressor parts only warranty Pyxis®, Omnicell® and AcuDose RX® compatible Exterior dimensions: 23 ¾"W x 24 ½"D x 33 ⅜"H Manual defrost 115V, 60 Hz, 1.5 Amps, 1/8 HP ENERGY STAR CERTIFIED Energy Consumption 0.97 kWh/day UL/C-UL listed Shipping weight: 131 lbs.</p>