

Product Name:	LRP-HCRFC9SS
Description:	<p>Plus Series Refrigerator &amp; Freezer Combination  9 Cu. Ft. Capacity (4.5 Ref./4.2 Freezer)  Two stainless steel exterior swing doors, right hinged*  Microprocessor temperature controller, refrigerator and freezer  Digital temperature display, refrigerator and freezer  Adjustable operating temperature range: Refrigerator: 1°C – 10°C;  Freezer:-15°C – -25°C  Audible and visual high and low temperature alarms, refrigerator and freezer  2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer  Remote alarm contacts, refrigerator and freezer  2 adjustable/1 fixed shelves refrigerator/2 fixed shelves freezer with guard rails on back in Refrigerator to prevent product from being pushed against back wall  Back wall evaporator cover, Refrigerator  Forced draft circulation Refrigerator  Keyed door locks  Magnetic door gasket for positive seal  Recessed door handles  Leveling legs  Stacking kit  LED interior lights  Probe access port (3/8"), 2 ea. (1 for refrigerator/1 for freezer)  Stainless steel exterior  HFC-free refrigerant (R600a), foam insulation and packaging  Two year parts and labor warranty, plus an additional three year compressor parts warranty  Exterior dimensions: 23 3/4" W x 24 1/2"D x 66 3/4"H  Manual defrost freezer / cycle defrost refrigerator  115V, 60 Hz, 1.74 Amps Refrigerator, 1.5 Amps Freezer, 1/8 HP Refrigerator and Freezer  UL/C-UL listed  Shipping weight: 266 lbs.  *Please note this unit requires 2 outlets.  *This product is available left hinged, product number LRP-HCRFC9SS-LH</p>