

Product Name:	LRP-KP-UCFS
Description:	Keyless Access Lock compatible with all LRP Plus Series Freestanding Undercounter Refrigerators & Freezers, and Compact Refrigerators Battery compartment is easily accessible from outside of the cabinet Includes 2-AAA batteries Convenient external battery installation/removal Low battery indicator Easy to clean keypad Must have 1 supervisor code and up to 19 user codes Factory install only Exterior Dimensions: 5.43" W x 1.13" D x 1.5" H Made in USA One year Manufacturer's Warranty