

Product Name:	LRP-MFB-30
Description:	<p>General Purpose Manual Defrost Laboratory Freezer  30 Cu. Ft. Capacity  One solid swing door, right hinged, not reversible  Mechanical thermostat  Adjustable operating temperature range: -15°-25°C  Four fixed shelves  Keyed door lock  Magnetic door gasket for positive seal  Stainless steel door handle  Leveling legs  Probe access port (3/8")  White textured cabinet and door  High density urethane foam cabinet and door insulation  CFC-free refrigerant (R404a), foam insulation and packaging  One year parts and labor warranty; plus additional four year  compressor parts warranty  Exterior dimensions: 35 ¾"W x 34 ¼"D x 72 ½"H  Manual defrost  115V, 60Hz, 6 Amps, 1/3 HP  UL/C-UL listed  Shipping weight: 365 lbs.</p>