

Product Name:	LRP-MFP-15-C
Description:	<p>General Purpose Series Manual Defrost Laboratory Chest Freezer</p> <p>15 Cu. Ft. capacity</p> <p>One solid door</p> <p>Microprocessor temperature controller:</p> <p>Adjustable operating temperature range: -15°C – -25°C</p> <p>Standard Digital Temperature Display and Alarm Module Sample simulated digital temperature display: .1C resolution, convertible to display temperatures in Celsius or Fahrenheit Audible and visual high/low temperature alarms</p> <p>Sensor error alarm</p> <p>Power failure alarm</p> <p>Resettable 1-point Min/Max History</p> <p>Remote alarm contacts</p> <p>LED Interior light</p> <p>Two removable baskets</p> <p>Swivel casters (2 ¾")</p> <p>Keyed door lock</p> <p>Magnetic door gasket for positive seal</p> <p>White cabinet and lid</p> <p>High density urethane foam cabinet and door insulation</p> <p>HFC-free refrigerant (R600a), foam insulation and packaging</p> <p>One-year parts and labor warranty, plus an additional four year compressor parts warranty</p> <p>Exterior dimensions: 56" W x 29 ½"D x 35 1/8"H with casters; 32 3/8" without casters (Height with hinged door open 90° is 62 ½")</p> <p>Manual defrost</p> <p>115V, 60 Hz, 5 Amps, 1/4 HP</p>

Shipping weight: 197 lbs.