

Product Name:	LRP-MFP-3030TS
Description:	<p>TempLog Plus Series Manual Defrost Laboratory Freezer, -30C  30 Cu. Ft. capacity  One solid swing door, right hinged, not reversible  Five inner "fast-freeze compartments" with doors for temperature maintenance and recovery  Microprocessor temperature controller:  Adjustable operating temperature range: -27°C – -33°C  Touch Screen Digital Temperature Display and Alarm Module  Password-protected control panel  Sample simulated digital temperature display; .1° Resolution, adjustable temperature display in Celsius or Fahrenheit  Active temperature scale bar display  Resettable 1-point min/max history  Digital 24-hour temperature chart  Temperature data logging  USB port for data transfer and software updates  Real-time operational display clock  4-20mA output  Electronic troubleshooting guide  Electronic QR Code to access Owner's Manual &amp; technical support contact Information  Adjustable alarm volume  Audible and visual high/low temperature alarms  Sensor error alarm  Power failure alarm  Door ajar alarm  Low battery alarm  Operation indicator  Alarm validation (audible, visual, and remote alarm)  Alarm mute with ring back  Remote alarm contacts  Four fixed shelves  Keyed door lock  Magnetic door gasket for positive seal  Stainless steel door handle  Caster platform (adds 4" to height)  One probe access port (3/8")  White textured cabinet and door  High density urethane foam cabinet and door insulation  CFC-free refrigerant (R404a), foam insulation and packaging  Three year parts and labor warranty, plus additional two year compressor parts warranty  Exterior dimensions: 35 3/4"W x 35"D x 82 1/2"H with caster platform; 78 1/2"H without caster platform  Manual defrost  115V, 60 Hz, 6 Amps, 1/3 HP  UL/C-UL listed  Shipping weight: 492 lbs.</p>