

Product Name:	PH-LRP-HC-16G
Description:	<p>Plus Series Pharmacy/Vaccine Glass Door Refrigerator  16 Cu. Ft. capacity  One swing glass door, self-closing (energy efficient, double pane vacuum insulated), right hinged, not reversible  Microprocessor Temperature Controller  Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C)  Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory. F/C switchable, field installable, and visual &amp; audible temp alarm  Standard Digital Temperature Display and Alarm Module with battery back-up to support temperature monitoring and alarms  Sample simulated digital temperature display; .1° Resolution  Adjustable temperature display in Celsius or Fahrenheit  Audible and visual high/low temperature alarms  Sensor error alarm  Power failure alarm  Resettable 1-point Min/Max History  Remote alarm contacts  Alarm validation (audible, visual, and remote alarm)  Alarm mute with ring back  2 probes; 1 in air and 1 in sample bottle  5 total shelves, 3 adjustable plus 2 additional  Forced draft air circulation  Keyed door lock  Vaccine storage power cord warning label  Magnetic door gasket for positive seal  Heavy duty door handle  Swivel Casters (3 ½”)  LED interior lights are safety shielded and switch controlled  Probe access port (3/8”)  Powder coated steel white interior and exterior  Hydrocarbon, natural refrigerant (R290)  Pharmacy refrigerator/freezer toolkit and temperature logs  Two year parts and labor warranty, plus an additional three year compressor parts warranty  Exterior dimensions: 25” W x 29 ¾” D x 79” H with casters; 75 ½” H without casters  Cycle defrost  115V, 60 Hz, 3 Amps, 1/5 HP  ENERGY STAR CERTIFIED  ETL listed  Shipping weight: 283 lbs</p>