



Exterior
PH-LRP-HC-RFC1020G


Front View



Side View



Top View



NOTE: Ventilation Requirements – Minimum 4" clearance on sides and back



115VAC/60Hz/1PH
NEMA-5-15P

Compartment	Cu. Ft	Defrost	Door	Int Door	Shelves	W"	D"	H"	Refrigerant	H.P.	Amps	Weight
Refrigerator Freezer	5.2 4	Cycle Manual	1 Glass 1 Solid	0 0	3 Adj. 2 Fixed	23 3/4	24	65	R600a	1/6 1/8	1.3 1.4	186 lbs.