



Exterior
PH-LRP-HC-RFC1040G


Front View



Side View



Top View



NOTE: Ventilation Requirements – Minimum 4" clearance on sides and back



115VAC/60Hz/1PH
NEMA-5-15P

Compartment	Cu. Ft	Defrost	Door	Int Door	Shelves	W"	D"	H"	Refrigerant	H.P.	Amps	Weight
Refrigerator Freezer	5.2 4	Cycle Manual	1 Glass 1 Solid	0 0	3 Adj. 3 Adj.	23 3/4	24	65	R600a R290	1/6 1/4	1.3 2.2	188 lbs.