

|               |  |
|---------------|--|
| Product Name: | PH-LRP-HC-RFC7   |
| Description:  | <p>Plus Series Pharmacy/Vaccine Refrigerator &amp; Freezer Combination<br/> 7 Cu. Ft. capacity (5.2 Ref./1.7 Freezer)<br/> Two exterior swing doors (glass door refrigerator/solid door freezer), right hinged, not reversible<br/> Microprocessor temperature controller, refrigerator and freezer<br/> Digital temperature display, refrigerator and freezer<br/> Adjustable operating temperature range: Refrigerator: 36°F – 46°F (2°C -8°C); Freezer: 5°F – -13°F (-15°C - -25°C)<br/> Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory, F/C switchable, field installable, and visual &amp; audible temp alarm<br/> Audible and visual high and low temperature alarms, refrigerator and freezer<br/> Remote alarm contacts, refrigerator and freezer<br/> 2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer<br/> 3 adjustable shelves refrigerator/2 adjustable shelves freezer, shelves have back rails to prevent product from being pushed against back wall<br/> Footed bottom shelf refrigerator, and freezer<br/> Vaccine storage power cord warning label<br/> Keyed door locks<br/> Magnetic door gasket for positive seal<br/> Integrated door handles<br/> Leveling legs<br/> Stacking kit<br/> Two probe access ports (3/8"); one, refrigerator/one, freezer<br/> White powder coated exterior<br/> HFC-free refrigerant (R600a), foam insulation and packaging<br/> Pharmacy refrigerator/freezer toolkit and temperature logs<br/> Two year parts and labor warranty, plus an additional three year compressor parts warranty<br/> Exterior dimensions: 23 5/8" W x 24 1/4"D x 53 1/8"H<br/> Manual defrost freezer / cycle defrost refrigerator<br/> 115V, 60 Hz, 1.3 Amps refrigerator, 1 Amps freezer, 1/6 HP refrigerator, 1/8 HP freezer<br/> UL/C-UL listed<br/> Shipping weight: 207 lbs.<br/> *Please note this unit requires 2 outlets.</p> |