

Product Name:	PH-LRP-HC-RFC7S
Description:	<p>Plus Series Pharmacy/Vaccine Refrigerator &amp; Freezer Combination  5.2 cu. ft. capacity (Refrigerator) / 1.7 cu. ft. capacity (Freezer)  Two exterior swing doors (solid door refrigerator/solid door freezer), right hinged, not reversible  Microprocessor temperature controller, refrigerator and freezer  Digital temperature display, refrigerator and freezer  Adjustable operating temperature range: Refrigerator: 36°F – 46°F (2°C -8°C); Freezer: 5°F – -13°F (-15°C - -25°C)  Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory. F/C switchable, field installable, and visual &amp; audible temp alarm  Audible and visual high and low temperature alarms, refrigerator and freezer  Remote alarm contacts, refrigerator and freezer  2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer  3 adjustable shelves refrigerator/2 adjustable shelves freezer, shelves have back rails to prevent product from being pushed against back wall  Footed bottom shelf refrigerator, and freezer  Forced draft air circulation, refrigerator  Vaccine storage power cord warning label  Keyed door locks  Magnetic door gasket for positive seal  Integrated door handles  Leveling legs  Stacking kit  Two probe access ports (3/8"); one, refrigerator/one, freezer  White powder coated exterior  HFC-free refrigerant (R600a), foam insulation and packaging  Pharmacy refrigerator/freezer toolkit and temperature logs  Two year parts and labor warranty, plus an additional three year compressor parts warranty  Exterior dimensions: 23 3/4"W x 24"D x 53 1/8"H  Manual defrost freezer / cycle defrost refrigerator  115V, 60 Hz, 1.3 Amps refrigerator, 1 Amp freezer, 1/6 HP refrigerator, 1/8 HP freezer  UL/C-UL listed  Shipping weight: 215 lbs.  *Please note this unit requires 2 outlets.</p>