

Product Name:	PH-LRP-HC-S72S
Description:	Basic Series Pharmacy/Vaccine Solid Door Refrigerator 72 Cu. Ft. capacity Three swing solid doors, self-closing, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts 2 probes; 1 in air and 1 in sample bottle Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C) Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory. F/C switchable, field installable, and visual & audible temp alarm Keyed door locks Heavy duty door handles Magnetic door gasket for positive seal Twelve adjustable shelves (adjustable in ½” increments) Vaccine storage power cord warning label Forced draft air circulation LED interior lights are safety shielded and switch controlled Swivel Casters (3 ½”) Probe access port (⅜”) White powder-coated steel exterior and interior Hydrocarbon, natural refrigerant (R290) Pharmacy refrigerator/freezer toolkit and temperature logs One year parts and labor warranty, plus an additional four year compressor parts warranty Exterior dimensions: 81” W x 34” D x 81 ¾” H with casters; 78 ¼” H without casters Cycle defrost 115V, 60 Hz, 7.5 Amps, 1/2 HP ETL listed Shipping weight: 600 lbs.