

Product Name:	LRP-HC-SSP-49
Description:	<p>Plus Series Stainless Steel Laboratory Refrigerator 49 Cu. Ft. Capacity Two stainless steel swing doors Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms Remote alarm contacts Adjustable operating temperature range: 1°C -10°C Six adjustable shelves Forced draft air circulation Keyed door locks Magnetic door gasket for positive seal Integrated door handles Swivel casters (adds 6 ½" to height) Interior light Probe access port (3/8") Stainless steel interior floor Stainless steel exterior High density urethane foam cabinet insulation Hydrocarbon, natural refrigerant (R290) Two year parts and labor warranty, plus an additional three year compressor parts warranty Exterior dimensions: 54"W x 35"D x 84"H with casters; 77 1/2" H without casters Cycle defrost 115V, 60 Hz, 4.5 Amps, 1/4 HP ETL listed Shipping weight: 450 lbs.</p>